




AALBORG UNIVERSITY

Testing of the Rehabilitation Robot

Additional material to group 363's semester report
ROB3 2020



Settling time is calculated by counting the measured steps from first correcting angle to first steady angle, and then dividing it with the frequency. This can be written as:

$$\text{settling time} = \frac{\text{steps}}{\text{frequency}}$$

This result is given in seconds.

Overshoot is calculated by taking the angle with the largest overshoot, dividing it by the steady state angle and multiplying it with 100. This is written as:

$$\text{overshoot} = \frac{\text{biggest overshoot sample}}{\text{steady state angle}} \cdot 100$$





This result is given in percentage.

Steady state error is calculated by subtracting the actual steady state angle by the desired steady state angle. This is written as:

$$\text{steady state error} = \text{actual angle} - \text{desired angle}$$

The result is given in degrees.

Below the results has been illustrated.



GUI

The interface will be tested by running a simulation of a normal rehabilitation, but without the physical aspect of the robot. One person will start the program, and every handling will be observed and described.

Handlings

- Clicking “Start”.
- Choses shoulder height between 3 options by clicking. Second option has been chosen.
- Choses distance to patience between 3 options by clicking. Third option has been chosen.
- Choses number of different positions (icebergs) between 3 options by clicking. First option has been chosen.
- “Done” is then pressed.
- A text box asks the user to place the positions within the map. Three positions are clicked on the screen.
- The game starts and can be played by the patient.
- The game is quit by clicking “Quit”

Overall clicks: 8-10 (depending on the number of positions)

Time between options exceeding 0.5 seconds: 0 (no lag)

Repetition

The repetition quality of the entire system will be tested by having the test person reach as many positions as possible within 1 minute. The positions (icebergs) will be placed with approximately 1/2 of the entire (diagonal) length of the map between them.

Results

Repetitions done: 5


To find the number of possible repetitions in 45 minutes the following calculations has been mad:

$5 \cdot 45 = 225$ repetitions


Precision

The precision of the system will be tested by having the test person follow a path close to the border of the map, at the same time the movement of the robot will be calculated by inverse kinematics. The position of the boat will be evaluated by Tracker 5.1.5 Video Analysis and Modeling Tool. The video will be analyzed at 18 FPS and 1080p. The two paths will then be compared.


Results



Horizontal is y-coordinates of the recorded movement and vertical is x-coordinates of the recorded movement, both in pixels.



Quick view of the video analysis tool



Output of inverse kinematics

Further information

Test done by: Group 363

Tested on: Benjamin Damsgaard (test person)

Date of test: 03/12/2020 - 04/12/2020

Place of testing: Fredrik Bajers Vej 7

Version: 2.0